

Rad-Sure™ ChromicVue  
15 Gy and 25 Gy  
Blood Irradiation Indicator  
Use Instructions  
Rev. 01 June 2019

When attached to blood products, Rad-Sure™ ChromicVue 15 Gy and 25 Gy indicators show whether the blood products have been irradiated. Before a blood product and its attached indicator are irradiated, the indicator reads “NOT IRRADIATED”. After the blood product and its attached indicator are irradiated, the word “NOT” in the indicator window is obscured and the indicator reads “IRRADIATED”.

Rad-Sure™ ChromicVue 15 Gy and 25 Gy indicators should be used only with irradiators having cesium-137 or cobalt-60 radiation sources, or other sources of equal or greater energy.

NOTE: Rad-Sure™ indicators only indicate that irradiation has occurred. They do not measure the dose from an irradiator. They must be used only with irradiators that are working properly and are properly calibrated. However, because of variation in lighting conditions, product batches and user eyesight, Rad-Sure™ indicators may appear fully opaque even after exposure to less than the minimum dose for the product. Therefore, Rad-Sure™ indicators should never be used to calibrate an irradiator.

1. Look at the top of the box for the arrow. Pull the clear strip to dispense one Rad-Sure™ indicator from the dispenser box.
2. Do not expose at or above 50°C (122°F)-Avoid prolonged exposure to elevated temperatures above 25°C.
3. Look for the word “NOT” in the window of each indicator. Do not use an indicator unless the word “NOT” is clearly visible (see figure 1a).
4. Observe the product expiration date on the indicator. If the product is beyond the expiration date, do not use the indicator.
5. Avoid touching the adhesive layer with rubber gloves. Apply the indicator firmly to an appropriate, clean, dry place on the unit of blood prior to irradiation.
6. Irradiate the unit of blood with the attached indicator in a gamma source irradiator at or above the recommended minimum dose.
7. After irradiation, check the indicator to verify that the window is black and obscures the word “NOT”. The indicator should now read “IRRADIATED” (see figure 1b). If the word “NOT” remains visible, notify your supervisor immediately.
8. Avoid prolonged exposure of the indicators to direct sunlight or ultraviolet light.
9. Store indicators at or below 25°C. Store away from all radiation sources, including gamma, rays, x-rays, electron devices, and microwaves.

Rad-Sure™ is a registered trademark owned by Ashland  
US Patent NO-7445880

Manufactured for:  
ASHLAND SPECIALTY INGREDIENTS G.P.  
1005 Rt. 202/206  
Bridgewater, NJ 08807  
USA

U.S. FDA Listing Number: A976736


For reorders within the U.S. email [AdvancedOrders@Ashland.com](mailto:AdvancedOrders@Ashland.com)

For service or information in the U.S. email: [AdvancedMaterials@Ashland.com](mailto:AdvancedMaterials@Ashland.com)


For service or information outside the U.S., contact your local distributor. For the name of a distributor in your area, contact: [AdvancedMaterials@Ashland.com](mailto:AdvancedMaterials@Ashland.com).

The information contained in this document is intended for use only by persons having technical skill and at their own discretion and risk. All information presented herein is believed to be accurate and reliable, but is not to be taken as a guarantee, an express warranty, or an implied warranty of merchantability or fitness for a particular purpose, or representation, express or implied, for which Ashland and its subsidiaries assume legal responsibility. Certain end uses of our products may be regulated pursuant to rules governing medical devices. A purchaser must make its own determination of our products' suitability for purchaser's use, whether alone or in conjunction with other products, the protection of the environment, and the health and safety of its employees and customers. It is also the purchaser's responsibility to comply with all applicable legal and regulatory requirements (including those governing medical devices). We make no warranty against infringement of any patents by reason of purchaser's use of any of our products.

**figure 1a**  
Rad-Sure™ ChromicVue™ 15 Gy and  
25 Gy indicators before irradiation







**figure 1b**  
Rad-Sure™ ChromicVue™ 15 Gy and  
25 Gy indicators after irradiation



## product list

product	product code
Rad-Sure™ ChromicVue™ 15 Gy blood irradiation indicators	943053
Rad-Sure™ ChromicVue™ 25 Gy blood irradiation indicators	943022

symbol/ símbolo/ symbole/ 符號/ Символ	symbol title/ nombre del símbolo/ explication du symbole/ symboltitel/ tytuł symbolu/ 符號說明/ Название символа
	Batch code/ Código de lote/ Code de lot/ Chargen-Code/ Kod partii/ 批次編號/ Код партии
	Use-by date/ Fecha de vencimiento/ Date limite d'utilisation/ Haltbarkeitsdatum/ Data ważności/ 到期/ Срок годности
	Manufacturer/ Fabricante/ Fabricant/ Hersteller/ Producent/ 製造商/ Производитель
	Do not re-use/ No reutilizar/ Ne pas réutiliser/ Darf nicht wiederverwendet werden/ Nie można go ponownie wykorzystywać/ 切勿重複使用/ Не использовать повторно