Klucel[™] Hydroxypropylcellulose in Hot-melt Extrusion Applications

Controlling drug solubility in immediate- and sustained-release tablets

Advantages of hot-melt extrusion (HME) over other processes include continuous processing efficiencies (high throughput, low cost, no solvents or water required), improved product uniformity, and enhanced solubility due to molecular dispersion of the active. HME can be used to make granules that are then tabletted or alternative dosage forms such as thin films (Figures 1 and 2).

In this application, molten thermoplastic polymers function as thermal binders during extrusion and upon cooling and solidification. Since most drugs are thermally sensitive, it is critical that the polymer can be processed at a relatively low temperature.

Klucel hydroxypropylcellulose (HPC) is the premier thermoplastic polymer for hot melt extrusion applications. This family of products is fully approved for pharmaceutical and food use and can be processed without plasticizers at significantly lower temperatures than other polymers (Table 1).

Low molecular weight grades such as Klucel EF and ELF can be used to enhance the solubility of low-soluble, Class II drugs. Klucel acts as a matrix in which drug molecules are immobilized and dispersed in a nanocrystralline or amorphous state.

Higher molecular weight grades (Klucel MF) can be used to make controlled-release tablets that are smaller, stronger, and less porous compared to those prepared via wet granulation. This allows a gel layer to be formed, retarding the dissolution of highly soluble drugs like Metformin (Figure 3).

Features and Benefits:

- Increased formulation stability nonionic polymer that does not require plasticizer
- Low temperature processing
- Enhanced solubility of low soluble actives
- High molecular weight grades for sustained release of highly soluble drugs
- Stronger, smaller, and less porous tablets compared to wet granulation
- Widest pharmaceutical and food regulatory approval compared to other thermoplastic polymers

Figure 1. Extruded strands of nifedipine, Klucel and MCC which have been milled and compressed into tablets.







Table 1. Typical glass transition temperatures of polymers commonly used in hot-melt extrusion.

Polymer	Tg*
PVP	165°C
Methylcellulose	150°C
HPMC (USP 2910 Type)	140°C
HPMCAS	120°C
EC	120°C
Copovidone (PVA-PVP)	115°C
Methacrylic acid copolymer E-type	125-160°C 50°C
НРС	0°C
PEO	-30°C



Global Headquarters

Ashland Inc. 50 East RiverCenter Blvd. P.O. Box 391 Covington, KY 41012-0391 U.S.A. Tel: +1 859 815 3333

Ashland Specialty Ingredients 8145 Blazer Drive Wilmington, DE 19808 U.S.A. Tel: +1 877 546 2782

Regional Centers

Europe — Switzerland Tel: +41 52 560 55 00 Fax: +41 52 560 55 99

India — Maharashtra Tel: +91 22 61489696 Fax: +91 22 61712514

Asia Pacific — Singapore Tel.: +65 6775 5366 Fax: +65 6776 9690

Latin America — Mexico Tel.: +5255 5211 0111 Fax: +5255 5212 0883

pharmaceutical@ashland.com

ashland.com/pharmaklucel

Registered trademark, Ashland or its subsidiaries, registered in various countrie ™ Trademark, Ashland or it subsidiaries, registered in various countries © 2010, 2013, Ashland © Conception Construction Construction Construction Construction Construction Construction Construction Construction © Construction Construction Construction Construction © Construction Construction Construction © Construc

All statements, information and data presented herein are believed to be accurate and reliable, but are not to be taken as a guarantee, an express warranty, or an implied warranty of merchantability or fitness for a particular purpose, or representation, express or implied, for which Ashland Inc. and its subsidiaries assume legal responsibility. The Ashland Research Center in Wilmington, Delaware is equipped with a Leistritz ZSE 18-mm twin-screw extruder that can be used for customer trials and new product development. We also have wide-ranging capabilities for thermal property analysis.

Please contact your local sales representative or visit us at www.aqualon.ashland.com/ pharmaceutical to learn more about Klucel HPC and its benefits in hot-melt extrusion.





* approximate: no plasticizers or other additives

Figure 3. Dissolution profile of extruded and wet granulated 75% Metformin / 25% Klucel MF tablets.



